

# HEALTH Monitor GMON



„INDICATE health risks – RATE professionally – REACT prophylactically“

## Register for scale TANITA Pro BCM Index card: measurement results BCM, Analysis

### Application

Connection to TANITA- Body composition analyzers

Electronic or manual data capture of all measurement results of the scale





### Evaluation of the results

Age – and specific evaluation of:

- BMI
- Body fat percentage and
- Body water level

Overview of the body composition and segment analysis (scale type MC-180)



### Report

Clear presentation of the measured values with evaluation of:

- BMI
- Muscle mass
- Body fat percentage and
- Body water percentage

Mass distribution (fat mass – Muscle mass)

Tabular form of the segment analysis divided into age- and gender specific reference ranges